






KC Haps  
Model 50.000




# Manual

**KC** Denmark


Research Equipment  
Limnology • Oceanography • Hydrobiology

	<p style="text-align: center;"><b>Manual for Haps</b></p>	<p style="text-align: center;"><b>Model no. 50.000</b></p>
	<div data-bbox="225 409 341 539" data-label="Image"> </div> <p>The Haps is very dangerous in unskilled hands and serious precautions must be taken to avoid accidents.</p> <p>KC-Denmark is not, and cannot be held, responsible for any damage(s) made to equipment or to operators who ignore safety precautions or because of misuse or wrong operation.</p>	<div data-bbox="954 315 1458 1014" data-label="Image"> </div>
	<p style="text-align: center;"><b>Preparing the Haps</b></p>	
<p>1</p>	<div data-bbox="225 1182 341 1312" data-label="Image"> </div> <p><b>IMPORTANT:</b> Be sure the automatic releaser is in lock position (lock grip in function). <u>Max. Operating depth is 200 m to avoid collapse of the releaser.</u></p> <p>You can extend the operating depth to full oceanographic depth (6000 m) by removal of the automatic releaser. When the Haps have been raised from the sea, it must remain on the winch and the sample tube must be removed <b><u>before lowering to the deck</u></b>. Otherwise the Haps will sink to the deck and prevent removal of the sample tube.</p>	<div data-bbox="959 1330 1453 1742" data-label="Image"> </div>



2	Place the bucket arm in locked position.	
3	Insert the sample tube.	
4	Mount the necessary amount of lead weights you need for penetrating the seabed.	



5	<p>When the Haps have been lowered from the sea, be sure the automatic releaser is in lock position. (Lock grip in function).</p> <p>If the releaser has been removed, the Haps must remain on the winch, see note #1.</p>	
6	<p>When removing the sample tube, a thin plate is placed under the sample.</p>	
7	<p>Remove the sample by slicing. The sample corer is inserted at top and each time you turn the handle 1 revolution, the sample will be released 5 mm, ready for slicing.</p>	

### Content of standard delivery

50.000	1	<p>Haps, frame supported bottom corer complete, consisting of:</p> <p>1 pc main rack, 4 pcs 10 kg lead weights (50.400). 4 pcs 5 kg lead weights. (50.200). 4 pcs 2 kg lead weights. (50.300). 1 pc non-return flap. (50.500) 1 pc sample tube, AISI 316 stainless steel. (50.600).</p>	108 kg		
--------	---	---	--------	---	--

### Other parts recommended

Order no.	Pcs.	Product	Net. weight		
50.100	10	Lead weight, 2 kg, for Haps, $\varnothing 150 \times \varnothing 36 \times 15$ mm. For mounting on top of the vertical $\varnothing 35$ mm rod.	2 kg each		
50.200	1	Lead weight, 5 kg, for Haps, $200 \times 100 \times 25$ mm. For mounting on the platform.	5 kg each		
50.300	20	Lead weight, 2 kg, for Haps, $\varnothing 150 \times \varnothing 13 \times 15$ mm. For mounting on the platform.	2 kg each		
50.400	4	10 kg lead weight for Haps, $540 \times 40 \times 40$ mm. For mounting on the bottom main frame.	10 kg each		
50.500	1	Non-return flap			
50.600	1	Sample tube, AISI 316 stainless steel. Length 315 mm, O.D./I.D. $\varnothing 140/\varnothing 136$ mm, sample area $145 \text{ cm}^2$ . Sample volume 4,5 litre.	3,5 kg		

50.700	1	Sample tube, polycarbonate. Length 315 mm, O.D./I.D. $\varnothing$ 133/ $\varnothing$ 127 mm, sample area 127 cm <sup>2</sup> . Sample volume 4,0 litre.	0,7 kg		
50.800	1	Sample ejection aggregate. Sample ejection aggregate  The sample corer is inserted at top and each time you turn the handle 1 revolution, the sample will be released 5 mm, ready for slicing.	10 kg		

Rev. March 23, 2005 – lkj

**KC** Denmark

Research Equipment  
Limnology • Oceanography • Hydrobiology

E-mail: [kc@kc-denmark.dk](mailto:kc@kc-denmark.dk) website: <http://www.kc-denmark.dk/>  
 Holmbladsvej 19 – DK-8600 Silkeborg – Denmark - Tel. +45 86 82 83 47 – Fax +45 86 82 49 50  
 Bank: Sydbank – S.W.I.F.T. SYBKDK22 IBAN DK5070460000104832 or  
 Bank: Danske Bank – S.W.I.F.T. DABADKKK IBAN6730003462779860  
 CVR no. DK 86 16 15 16